



**CLINICAL
MICROBIOLOGY**

Index.

Sample processing.	4
<hr/>	
Transportation	4
Pretreatment	4
Clinical diagnostic.	5
<hr/>	
General Analysis	5
Anaerobic Analysis	6
Mycobacteria Analysis	6
Specific Analysis	6
Digestive tract	
Urinary tract	
Respiratory tract	
Genital tract	
Blood Cultures	
Mycoplasmas	
CondaChrome®	
Supplements	8
Antimicrobial susceptibility.	9
<hr/>	
Culture media	9
CondaChrome®	10
MIC	10
Monitoring.	22
<hr/>	
Environmental	22
Medical devices	22
Stains and dyes.	23
<hr/>	
Other reagents.	24
<hr/>	
Ingredients	24
Agaroses	24
Strain Preservation	24
Antibiotics	24

Clinical microbiology and public health


Clinical microbiology is a tool used to isolate and identify microorganisms, generally of the pathogenic kind, with the aim of providing a diagnosis and suggesting a treatment to patients.

The close contact between clinical microbiology and public health is fundamental to:

- The **prevention and management** of infectious diseases
- The **monitoring** of local and global epidemics
- The development and adoption of **new technologies and therapeutic methods**


Globally, out of a 100 hospitalized patients, between 7 and 10 will suffer at least one nosocomial infection. These infections are associated with high morbidity and mortality, and in the U.S. alone are responsible for more than 90,000 deaths per year, according to the Center for Disease Control (CDC). This is very similar in Europe, where 6 nosocomial infections are the main cause of death, according to the ECDC.

What are the main culprits of nosocomial infections?



When are nosocomial infections acquired?

23% of Healthcare-associated infections (HAIs) are detected on **admission**, of which 54% are associated with a previous stay in the **same hospital**.





Sample Processing.

The work of the clinical microbiology laboratory involves paying special attention even prior to the sample analysis, which begins by ensuring an adequate sample collection using appropriate containers and culture media for collection.

Likewise, in order to guarantee the preservation and integrity of the samples during the transportation, it is necessary to use media specifically developed for this purpose. And due to the nature of some clinical samples, even a pre-treatment is necessary, consisting on the use of certain culture media and different separation or homogenization techniques.

CATEGORY	CAT.	PRODUCT	PACK SIZE
Pretreatment	1000	Anaerobic Agar	500 g
Pretreatment	1052	Macconkey Agar	500 g
Pretreatment	1090	Sabouraud Dextrose Agar + Chloramphenicol	500 g
Pretreatment	1106	GC Agar Base	500 g
Pretreatment	1107	Bordet-Gengou Agar Base	500 g
Pretreatment	1108	Blood Agar Base	500 g
Pretreatment	1117	Campylobacter Medium Base	500 g
Pretreatment	1152	Columbia CNA Agar Base	500 g
Pretreatment	1400	Brain Heart Infusion Broth (BHI Broth)	500 g
Pretreatment	1508	Thioglycollate Fluid Medium	500 g
Pretreatment	4220	Roiron Medium	20 Tubes
Pretreatment	5015	Reagent N° 2	10 x 100 ml
Pretreatment	5017	B Reagent	10 x 100 ml
Pretreatment	6006	Reagent N°1	10 capsules
Transportation	1236	Todd Hewitt Broth	500 g
Transportation	1518	Stuart Transport Medium	500 g
Transportation	1529	Cary-Blair Medium	500 g
Transportation	1530	Amies Transport Medium without Charcoal	500 g
Transportation	1535	Amies Transport with Charcoal	500 g



Clinical Diagnostic.

Condalab offers a wide range of culture media for routine clinical diagnosis, CE/IVD Marking, as well as specific products for certain nosocomial pathogens, and a wide range of CondaChrome® media for rapid microbial identification.

Our products meet the sterility conditions that are so important in the diagnosis of infectious diseases. Their easy use, combined with an optimal recovery and differentiation of pathogens make our culture media and reagents the best ally for daily clinical practice.

General Analysis

CAT.	PRODUCT	PACK SIZE
1024	Sabouraud Dextrose Agar	500 g
1028	DNAse Test Agar (Deoxyribonuclease Activity)	500 g
1031	Bile Esculin Agar	500 g
1039	Eosin Methylene Blue Agar (EMB)	500 g
1048	Brain Heart Infusion Agar (BHI Agar)	500 g
1050	Levine Agar (EMB)	500 g
1052	Macconkey Agar	500 g
1058	Mueller Hinton Agar	500 g
1090	Sabouraud Dextrose Agar + Chloramphenicol	500 g
1104	Columbia Agar Base	500 g
1106	GC Agar Base	500 g
1108	Blood Agar Base	500 g
1152	Columbia CNA Agar Base	500 g
1224	Trypticasein Soy Broth (TSB)	500 g
1400	Brain Heart Infusion Broth (BHI Broth)	500 g
1500	OF Basal Medium (Hugh And Leifson)	500 g
1508	Thioglycollate Fluid Medium	500 g
1551	Luria Broth (Miller's LB Broth)	500 g
2059	Sabouraud Dextrose Agar + Chloramphenicol + Gentamicin	500 g
2180	Urea Agar Base (Christensen)	500 g

Anaerobic Analysis

CAT.	PRODUCT	PACK SIZE
1000	Anaerobic Agar	500 g
1012	Brucella Agar	500 g
1046	Triple Sugar Iron Agar (TSI)	500 g
1066	Schaedler Agar	500 g
1447	Clostridium Difficile Agar Base	500 g

Mycobacteria Analysis

CAT.	PRODUCT	PACK SIZE
1037	Macconkey Agar w/o Crystal Violet	500 g
1226	Urea Broth	500 g
1439	Lowenstein Jensen Medium Base w/o M.Green	500 g
2042	Middlebrook Broth	500 g
2043	Middlebrook Agar	500 g
4100	Lowenstein Jensen Medium	20 Tubes
4120	Lowenstein Jensen Medium with Pyruvate	20 Tubes
4125	Lowenstein Jensen Medium (Long Tube)	20 Tubes
4126	Lowenstein Jensen Medium with Pyruvate (Long Tube)	20 Tubes
4127	Coletsos Medium	20 Tubes

Specific Analysis

CATEGORY	CAT.	PRODUCT	PACK SIZE
Conda Chrome ®	1122	Salmonella Chromogenic Agar	500 g
Conda Chrome ®	1382	Candida Chromogenic Agar	500 g
Conda Chrome ®	1493	Pseudomonas Chromogenic Agar	500 g
Conda Chrome ®	1588	E. coli O157:H7 Cromogenic Agar Base	500 g
Conda Chrome ®	2054	Vibrio Chromogenic Agar	500 g
Conda Chrome ®	2076	Staphylococcus Chromogenic Agar	500 g
Conda Chrome ®	2119	Klebsiella Chromogenic Agar Base	500 g
Blood Cultures	3001	Trypticasein Soy Broth with SPS, CO2 & Vacuum	10 x 50 ml
Blood Cultures	3004	Brain Heart Infusion Broth with SPS, CO2 & Vacuum	10 x 50 ml
Blood Cultures	3005	Brain Heart Infusion Broth with SPS, CO2 & Vacuum	8 x 20 ml
Blood Cultures	3008	Trypticasein Soy Broth with SPS, CO2 & Vacuum	8 x 20 ml
Blood Cultures	3105	Thioglycollate Broth with SPS, CO2 & Vacuum	10 x 50 ml
Blood Cultures	3107	Schaedler Broth with SPS, Vaccum and CO2	10 x 50 ml
Mycoplasmas	1140	PPLO Agar Base w/o Crystal Violet	500 g
Mycoplasmas	1262	PPLO Broth Base w/o Crystal Violet	500 g
Digestive tract	1011	Bismuth Sulfite Agar (Wilson Blair)	500 g
Digestive tract	1030	Hektoen Enteric Agar	500 g
Digestive tract	1044	Lysine Iron Agar	500 g
Digestive tract	1064	Salmonella Shigella Agar (SS Agar)	500 g
Digestive tract	1074	TCBS Agar	500 g

CATEGORY	CAT.	PRODUCT	PACK SIZE
Digestive tract	1078	Brilliant Green Agar	500 g
Digestive tract	1080	XLD Agar (Xylose Lysine Desoxycholate Agar)	500 g
Digestive tract	1099	Macconkey Agar with Sorbitol (CT-SMAC)	500 g
Digestive tract	1110	Urea Agar Base (Christensen)	500 g
Digestive tract	1114	Tetrathionate Broth Base	500 g
Digestive tract	1126	Yersinia Selective Agar Base (CIN)	500 g
Digestive tract	1129	Campylobacter Agar Base Blood Free (CCDA)	500 g
Digestive tract	1130	Muller Kauffmann Tetrathionate Broth Base	500 g
Digestive tract	1131	Campylobacter Agar Base (Preston)	500 g
Digestive tract	1186	Modified Salmonella Shigella Agar	500 g
Digestive tract	1222	Sodium Selenite Broth	500 g
Digestive tract	1248	GN Enrichment Broth (Hajna)	500 g
Digestive tract	1298	Sorbitol Peptone Broth And Bile Salts (PSB)	500 g
Digestive tract	1314	Nutrient Agar Dev Regulations	500 g
Digestive tract	1360	Salmonella Shigella Agar with Sodium Deoxycholate and Calcium Chloride (SSDC)	500 g
Digestive tract	1364	Kligler Iron Agar	500 g
Digestive tract	1370	Agar Base Aeromonas (RYAN)	500 g
Digestive tract	1401	Buffered Peptone Water	500 g
Digestive tract	1407	Alkaline Peptone Water	500 g
Digestive tract	1414	Rappaport Vassiliadis Broth	500 g
Digestive tract	1503	Wilkins Chalgren Medium	500 g
Digestive tract	1568	Wilkins Chalgren Medium II	500 g
Digestive tract	2024	Agar Test pH 8,0 for Inhibition Test	500 g
Digestive tract	2066	Helicobacter Agar Modified	500 g
Digestive tract	2097	Inositol Brilliant Green Bile Agar (Plesiomonas Differential Agar)	500 g
Genital tract	1060	Nutrient Agar	500 g
Genital tract	2036	Granada Agar Base	500 g
Genital tract	2099	Streptococcus Selective Agar	500 g
Respiratory tract	1022	Potato Dextrose Agar	500 g
Respiratory tract	1062	Mannitol Salt Agar (MSA) (Chapman Medium)	500 g
Respiratory tract	1072	Mycobiotic Agar (Fungal Selective Agar)	500 g
Respiratory tract	1088	Sabouraud Dextrose Agar + Cycloheximide	500 g
Respiratory tract	1089	Sabouraud Dextrose Agar + Chloramphenicol + Cycloheximide	500 g
Respiratory tract	1102	Cetrimide Agar Base	500 g
Respiratory tract	1134	Sabouraud Dextrose Agar + Chloramphenicol	500 g
Respiratory tract	1311	Legionella BCYE Agar Base	500 g
Respiratory tract	1317	Tinsdale Agar Base	500 g
Respiratory tract	1347	Burkholderia Cepacia Agar Base	500 g
Respiratory tract	1546	YPD Agar	500 g
Respiratory tract	1557	Legionella Agar Base + ACES + Alpha Ketoglutaric + Ferric pyrophosphate	500 g
Respiratory tract	2087	Loeffler Serum Medium Base	500 g
Urinary tract	1016	CLED Agar (Cystine Lactose Electrolyte Deficient)	500 g
Urinary tract	1068	Trypticasein Soy Agar (TSA)	500 g
Urinary tract	1424	Urinary Tract Infections Chromogenic Agar (UTIC)	500 g

Supplements

CAT.	PRODUCT	PACK SIZE
5100	Urea 40% Solution	100 ml
5152	Egg Yolk Emulsion	100 ml
6011	Polyenrichment Supplement	2x5 vials
6012	LCAT Supplement	10 vials
6013	VCN Supplement	10 vials
6014	VCAT Supplement	10 vials
6019	Preston Campylobacter Supplement	10 vials
6022	Legionella BCYE Supplement	10 vials
6025	Legionella GVPC Supplement	10 vials
6026	VCNT Supplement	10 vials
6032	Burkholderia cepacia Supplement	10 vials
6033	Yersinia Selective Supplement (CIN)	10 vials
6037	OADC Enrichment Supplement	10 vials
6038	ADC Enrichment Supplement	10 vials
6052	Ampicillin Supplement (Aeromonas)	10 vials
6053	CCDA Supplement (Campylobacter Blood free)	10 vials
6061	Clostridium Difficile Supplement	10 vials



Antimicrobial Sensitivity.

The analysis of microorganism sensitivity to different antimicrobial agents is a key point in any clinical microbiology laboratory. This analysis is carried out by the use of an antibiogram, the objective of which is to evaluate in a laboratory the response of a microorganism to a battery of antibiotics or antifungals.

Condalab offers classic antibiotic discs and strips with concentration gradients under EUCAST and CLSI control, as well as a specific line of CondaChrome® media for the detection of resistant strains. The monitoring of antimicrobial sensitivity and resistance is extremely important in the clinical field to provide the most effective treatment.

CATEGORY	CAT.	PRODUCT	PACK SIZE
Culture media	1015	Czapek-Dox Modified Agar	500 g
Culture media	1022	Potato Dextrose Agar	500 g
Culture media	1024	Sabouraud Dextrose Agar	500 g
Culture media	1046	Triple Sugar Iron Agar (TSI)	500 g
Culture media	1048	Brain Heart Infusion Agar (BHI Agar)	500 g
Culture media	1052	Macconkey Agar	500 g
Culture media	1058	Mueller Hinton Agar	500 g
Culture media	1062	Mannitol Salt Agar (MSA) (Chapman Medium)	500 g
Culture media	1088	Sabouraud Dextrose Agar + Cycloheximide	500 g
Culture media	1089	Sabouraud Dextrose Agar + Chloramphenicol + Cycloheximide	500 g
Culture media	1090	Sabouraud Dextrose Agar + Chloramphenicol	500 g
Culture media	1108	Blood Agar Base	500 g
Culture media	1134	Sabouraud Dextrose Agar + Chloramphenicol	500 g
Culture media	1152	Columbia CNA Agar Base	500 g
Culture media	1157	Dermatophytes Agar Base	500 g
Culture media	1160	Rose Bengal Agar + Chloramphenicol + Dichloran (DRBC Agar)	500 g
Culture media	1161	Dichloran Glycerol Agar (DG 18)	500 g
Culture media	1164	Corn Meal Agar	500 g
Culture media	1214	Mueller Hinton Broth	500 g
Culture media	1224	Trypticasein Soy Broth (TSB)	500 g
Culture media	1301	Chloramphenicol Agar (YGC Agar)	500 g
Culture media	1372	Bile Esculin Azide Agar	500 g
Culture media	1400	Brain Heart Infusion Broth (BHI Broth)	500 g
Culture media	1500	OF Basal Medium (Hugh And Leifson)	500 g
Culture media	1527	OGA Medium (Oxytetracycline Glucose Agar Base) (OGYE)	500 g

CATEGORY	CAT.	PRODUCT	PACK SIZE
Culture media	1545	Yeast Nitrogen Base w/o Amino Acids	500 g
Culture media	1546	YPD Agar	500 g
Culture media	1547	YPD Broth	500 g
Culture media	2042	Middlebrook Broth	500 g
Culture media	2116	Motility Test Medium	500 g
Culture media	2180	Urea Agar Base (Christensen)	500 g
Conda Chrome ®	1382	Candida Chromogenic Agar	500 g
Conda Chrome ®	1423	MRSA Chromogenic Agar Base	500 g
Conda Chrome ®	1498	MRSA Chromogenic Modified Agar Base	500 g
Conda Chrome ®	2062	ESBL Chromogenic Agar	500 g
Conda Chrome ®	2063	KPC Chromogenic Medium	500 g
Conda Chrome ®	2076	Staphylococcus Chromogenic Agar	500 g
Conda Chrome ®	2077	Vancomycin-Resistant Enterococcus (VRE) Chromogenic Agar	500 g
MIC	7000	Amikacin 30 µg	5x50 Discs
MIC	7001	Amoxicillin 25 µg	5x50 Discs
MIC	7002	Amoxicillin / Clavulanic acid 30 µg	5x50 Discs
MIC	7003	Amphotericin B 100 U µg	5x50 Discs
MIC	7004	Ampicillin 10 µg	5x50 Discs
MIC	7005	Ampicillin / Sulbactam 20 µg	5x50 Discs
MIC	7006	Azithromycin 15 µg	5x50 Discs
MIC	7007	Aztreonam 30 µg	5x50 Discs
MIC	7008	Bacitracin 0.04 µg	5x50 Discs
MIC	7009	Cefaclor 30 µg	5x50 Discs
MIC	7010	Cefamandole 30 µg	5x50 Discs
MIC	7011	Cefazolin 30 µg	5x50 Discs
MIC	7012	Cefepime 30 µg	5x50 Discs
MIC	7013	Cefepime / Clavulanic acid 40 µg	5x50 Discs
MIC	7014	Cefixime 30 µg	5x50 Discs
MIC	7015	Cefixime 5 µg	5x50 Discs
MIC	7016	Cefoperazone 75 µg	5x50 Discs
MIC	7017	Cefotaxime 30 µg	5x50 Discs
MIC	7018	Cefotaxime 5 µg	5x50 Discs
MIC	7019	Cefotaxime / Boronic acid 430 µg	5x50 Discs
MIC	7020	Cefotaxime / Clavulanic acid 40 µg	5x50 Discs
MIC	7021	Cefotetan 30 µg	5x50 Discs
MIC	7022	Cefoxitin 30 µg	5x50 Discs
MIC	7023	Cefpirome 30 µg	5x50 Discs
MIC	7024	Cefpodoxime 10 µg	5x50 Discs
MIC	7025	Cefsulodin 30 µg	5x50 Discs
MIC	7026	Ceftazidime 30 µg	5x50 Discs
MIC	7027	Ceftazidime / Clavulanic acid 40 µg	5x50 Discs
MIC	7028	Ceftizoxime 30 µg	5x50 Discs
MIC	7029	Ceftriaxone 30 µg	5x50 Discs
MIC	7030	Cefuroxime 30 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7031	Cephalexin 30 µg	5x50 Discs
MIC	7032	Cephalothin 30 µg	5x50 Discs
MIC	7033	Chloramphenicol 30 µg	5x50 Discs
MIC	7034	Ciprofloxacin 5 µg	5x50 Discs
MIC	7035	Clarithromycin 15 µg	5x50 Discs
MIC	7036	Amikacin 10 µg	5x50 Discs
MIC	7037	Aminosidin 60 µg	5x50 Discs
MIC	7038	Amoxicillin 10 µg	5x50 Discs
MIC	7039	Amoxicillin 2 µg	5x50 Discs
MIC	7040	Clindamycin 2 µg	5x50 Discs
MIC	7041	Clotrimazole 10 µg	5x50 Discs
MIC	7042	Trimethoprim / Sulphamethoxazol 25 µg	5x50 Discs
MIC	7043	Colistin 10 µg	5x50 Discs
MIC	7044	Doxyxyciclina 30 ug	5x50 Discs
MIC	7045	Econazole 50 µg	5x50 Discs
MIC	7046	Trimethoprim 25 µg	5x50 Discs
MIC	7047	Ertapenem 10 µg	5x50 Discs
MIC	7048	Erythromycin 15 µg	5x50 Discs
MIC	7049	Fosfomycin 200 µg	5x50 Discs
MIC	7050	Fusidic acid 10 µg	5x50 Discs
MIC	7051	Gentamicin 10 µg	5x50 Discs
MIC	7052	Imipenem 10 µg	5x50 Discs
MIC	7053	Kanamycin 30 µg	5x50 Discs
MIC	7054	Kanamycin 1000 µg	5x50 Discs
MIC	7055	Ketoconazole 50 µg	5x50 Discs
MIC	7056	Levofloxacin 5 µg	5x50 Discs
MIC	7057	Lincomycin 15 µg	5x50 Discs
MIC	7058	Mecillinam 10 µg	5x50 Discs
MIC	7059	Meropenem 10 µg	5x50 Discs
MIC	7060	Miconazole 10 µg	5x50 Discs
MIC	7061	Minocycline 30 µg	5x50 Discs
MIC	7062	Moxifloxacin 5 µg	5x50 Discs
MIC	7063	Nalidixic acid 30 µg	5x50 Discs
MIC	7064	Neomycin 30 µg	5x50 Discs
MIC	7065	Netilmicin 30 µg	5x50 Discs
MIC	7066	Nitrofurantoin 300 µg	5x50 Discs
MIC	7067	Norfloxacin 10 µg	5x50 Discs
MIC	7068	Novobiocin 5 µg	5x50 Discs
MIC	7069	Nystatin 100 U µg	5x50 Discs
MIC	7070	Ofloxacin 10 µg	5x50 Discs
MIC	7071	Ofloxacin 5 µg	5x50 Discs
MIC	7072	ONPG	5x50 Discs
MIC	7073	Optochin 5 µg	5x50 Discs
MIC	7074	Oxacillin 1 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7075	Oxidase Strip	50 strips
MIC	7076	Pefloxacin 5 µg	5x50 Discs
MIC	7077	Penicillin G 10 U	5x50 Discs
MIC	7078	Piperacillin 100 µg	5x50 Discs
MIC	7079	Piperacillin / Tazobactam 110 µg	5x50 Discs
MIC	7080	Polymixin B 300 U µg	5x50 Discs
MIC	7081	Rifampin (Rifampicin) 5 µg	5x50 Discs
MIC	7082	Spectinomycin 100 µg	5x50 Discs
MIC	7083	Streptomycin 10 µg	5x50 Discs
MIC	7084	Streptomycin 300 µg	5x50 Discs
MIC	7085	Amoxicillin 30 µg	5x50 Discs
MIC	7086	Amoxicillin / Clavulanic acid 3 µg	5x50 Discs
MIC	7087	Amphotericin B 20 U	5x50 Discs
MIC	7088	Teicoplanin 30 µg	5x50 Discs
MIC	7090	Tetracycline 25 µg	5x50 Discs
MIC	7091	Tetracycline 30 µg	5x50 Discs
MIC	7092	Ticarillin 75 µg	5x50 Discs
MIC	7093	Tigecycline 15 µg	5x50 Discs
MIC	7094	Tobramycin 10 µg	5x50 Discs
MIC	7095	Ampicillin 30 µg	5x50 Discs
MIC	7096	Azlocillin 75 µg	5x50 Discs
MIC	7097	Carbenicillin 100 µg	5x50 Discs
MIC	7098	Cefadroxil 30 µg	5x50 Discs
MIC	7099	Cefdinir 5 µg	5x50 Discs
MIC	7100	Cefonicid 30 µg	5x50 Discs
MIC	7101	Cefotaxime 10 µg	5x50 Discs
MIC	7102	Cefprozil 30 µg	5x50 Discs
MIC	7103	Ceftazidime 10 µg	5x50 Discs
MIC	7104	Ceftibuten 30 µg	5x50 Discs
MIC	7106	Ciprofloxacin 1 µg	5x50 Discs
MIC	7107	Clindamycin 10 µg	5x50 Discs
MIC	7108	Cloxacillin 1 µg	5x50 Discs
MIC	7110	Doripenem 10 µg	5x50 Discs
MIC	7111	Enoxacin 10 µg	5x50 Discs
MIC	7112	Fluconazole 10 µg	5x50 Discs
MIC	7113	Flucytosine 1 µg	5x50 Discs
MIC	7114	Fosfomycin 50 µg	5x50 Discs
MIC	7115	Gatifloxacin 5 µg	5x50 Discs
MIC	7116	Gemifloxacin 5 µg	5x50 Discs
MIC	7117	Gentamicin 120 µg	5x50 Discs
MIC	7118	Gentamicin 30 µg	5x50 Discs
MIC	7119	Gentamicin 500 µg	5x50 Discs
MIC	7120	Itraconazole 10 µg	5x50 Discs
MIC	7121	Linezolid 30 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7122	Linezolid 10 µg	5x50 Discs
MIC	7123	Lomefloxacin 10 µg	5x50 Discs
MIC	7124	Metronidazole 5 µg	5x50 Discs
MIC	7125	Mezlocillin 75 µg	5x50 Discs
MIC	7126	Mezlocillin 30 µg	5x50 Discs
MIC	7127	Moxalactam 30 µg	5x50 Discs
MIC	7128	Mupirocin 5 µg	5x50 Discs
MIC	7129	Mupirocin 200 µg	5x50 Discs
MIC	7130	Nafcillin 1 µg	5x50 Discs
MIC	7131	Nitrofurantoin 100 µg	5x50 Discs
MIC	7132	Nitrofurazone 100 µg	5x50 Discs
MIC	7133	Norfloxacin 5 µg	5x50 Discs
MIC	7134	Penicillin G 1 U	5x50 Discs
MIC	7135	Piperacillin 30 µg	5x50 Discs
MIC	7136	Piperacillin / Tazobactam 36 µg	5x50 Discs
MIC	7137	Ampicillin 2 µg	5x50 Discs
MIC	7138	Sparfloxacin 5 µg	5x50 Discs
MIC	7139	Spiramycin 10 µg	5x50 Discs
MIC	7140	Streptomycin 500 µg	5x50 Discs
MIC	7141	Sulfisoxazole 250 µg	5x50 Discs
MIC	7142	Ampicillin 25 µg	5x50 Discs
MIC	7144	Ticarillin / Clavulanic acid 85 µg	5x50 Discs
MIC	7145	Vancomycin 30 µg	5x50 Discs
MIC	7146	Oxidase Test Discs	5x50 Discs
MIC	7147	Dirithromycin 15 µg	5x50 Discs
MIC	7149	Ampicillin/cloxacillin 30 µg	5x50 Discs
MIC	7151	Ampicillin / Dicloxacillin 11 µg	5x50 Discs
MIC	7152	Ampicillin / Sulbactam 30 µg	5x50 Discs
MIC	7153	Apramycin 15 µg	5x50 Discs
MIC	7154	Apramycin 30 µg 5x50	5x50 Discs
MIC	7155	Azithromycin 30 µg	5x50 Discs
MIC	7156	Aztreonam 10 µg	5x50 Discs
MIC	7157	Bacitracin 10 U	5x50 Discs
MIC	7158	Dispenser (6 Magazine)	1 Dispenser
MIC	7159	X + V factor disc	5x50 Discs
MIC	7160	X Factor Disc	5x50 Discs
MIC	7161	Boronic acid 250 µg	5x50 Discs
MIC	7162	Boronic acid 300 µg	5x50 Discs
MIC	7163	Boronic acid 400 µg	5x50 Discs
MIC	7164	Carbenicillin 25 µg	5x50 Discs
MIC	7165	Cefaclor 10 µg	5x50 Discs
MIC	7166	Cefdinir 30 µg	5x50 Discs
MIC	7167	Cefditoren 5 µg	5x50 Discs
MIC	7168	Cefixime 10 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7169	Cefodizime 30 µg	5x50 Discs
MIC	7170	Cefoperazone 30 µg	5x50 Discs
MIC	7171	Cefoperazone / Sulbactam 105 µg	5x50 Discs
MIC	7172	Prulifloxacin 5 µg	5x50 Discs
MIC	7173	Vancomycin 120 µg	5x50 Discs
MIC	7174	Cefoxitin 10 µg	5x50 Discs
MIC	7175	Cefoxitin / Boronic acid 430 µg	5x50 Discs
MIC	7176	Cefpodoxime 30 ug	5x50 Discs
MIC	7177	Cefpodoxime / Clavulanic acid 11 µg	5x50 Discs
MIC	7178	Cefpodoxime / Clavulanic acid 40 µg	5x50 Discs
MIC	7179	Cefquinome 10 µg	5x50 Discs
MIC	7180	Cefquinome 30 µg	5x50 Discs
MIC	7182	Ceftazidime / Boronic acid 430 µg	5x50 Discs
MIC	7183	Ceftazidime / Tazobactam 40 µg	5x50 Discs
MIC	7184	Ceftiofur 30 µg	5x50 Discs
MIC	7185	Ceftriaxone 10 µg	5x50 Discs
MIC	7186	Ceftriaxone 5 µg	5x50 Discs
MIC	7187	Ceftriaxone / Clavulanic acid 40 ug	5x50 Discs
MIC	7188	Cefuroxime 5 µg	5x50 Discs
MIC	7189	Cefuroxime axetil 10 µg	5x50 Discs
MIC	7190	Cefuroxime axetil 30 µg	5x50 Discs
MIC	7191	Cephalonium 30 µg	5x50 Discs
MIC	7192	Cephradine 30 µg	5x50 Discs
MIC	7193	Chloramphenicol 10 µg	5x50 Discs
MIC	7194	Chloramphenicol 5 µg	5x50 Discs
MIC	7195	Chloramphenicol 50 µg	5x50 Discs
MIC	7196	Chlortetracycline 30 µg	5x50 Discs
MIC	7197	Ciclopirox olamine 50 µg	5x50 Discs
MIC	7198	Imipenem / Cilastatin 20 µg	5x50 Discs
MIC	7199	Ciprofloxacin 10 µg	5x50 Discs
MIC	7200	Ciprofloxacin 30 µg	5x50 Discs
MIC	7201	Ciprofloxacin 50 µg	5x50 Discs
MIC	7202	Clarithromycin 2 µg 5 µg	5x50 Discs
MIC	7203	Clarithromycin 5 µg	5x50 Discs
MIC	7204	Clinafloxacin 5 ug	5x50 Discs
MIC	7205	Clindamycin 15 µg	5x50 Discs
MIC	7206	Clotrimazole 300 µg	5x50 Discs
MIC	7207	Clotrimazole 50 µg	5x50 Discs
MIC	7208	Cloxacillin 10 µg	5x50 Discs
MIC	7209	Cloxacillin 30 µg	5x50 Discs
MIC	7210	Cloxacillin 5 µg	5x50 Discs
MIC	7211	Cloxacillin 5 µg 5x50 Discs	5x50 Discs
MIC	7212	Colistin 100 µg	5x50 Discs
MIC	7213	Colistin 25 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7214	Colistin 50 µg	5x50 Discs
MIC	7215	Danofloxacin 5 µg	5x50 Discs
MIC	7216	Daptomycin 30 µg	5x50 Discs
MIC	7217	Dibekacin 10 ug	5x50 Discs
MIC	7218	Difloxacin 10 ug	5x50 Discs
MIC	7219	Dispenser (8 Magazine)	1 Dispenser
MIC	7220	Doxycycline 10 µg	5x50 Discs
MIC	7221	Econazole 10 µg	5x50 Discs
MIC	7222	EDTA 750 µg	5x50 Discs
MIC	7223	EDTA 930 µg	5x50 Discs
MIC	7224	Ertapenem / EDTA 292 µg	5x50 Discs
MIC	7225	Ertapenem/EDTA 750 µg 5 µg	5x50 Discs
MIC	7226	Imipenem / EDTA 292 µg	5x50 Discs
MIC	7227	Imipenem/ EDTA 10 µg	5x50 Discs
MIC	7228	Imipenem / EDTA 750 µg	5x50 Discs
MIC	7229	Meropenem / EDTA 292 µg	5x50 Discs
MIC	7230	Meropenem/EDTA 750 µg	5x50 Discs
MIC	7231	EDTA / SMA 750 µg	5x50 Discs
MIC	7232	Enrofloxacin 10 µg	5x50 Discs
MIC	7233	Enrofloxacin 5 µg	5x50 Discs
MIC	7234	Erythromycin 10 µg	5x50 Discs
MIC	7235	Erythromycin 30 µg	5x50 Discs
MIC	7236	Erythromycin 5 µg	5x50 Discs
MIC	7237	Fleroxacin 5 µg	5x50 Discs
MIC	7238	Florfenicol 30 µg	5x50 Discs
MIC	7239	Flucloxacillin 5 µg	5x50 Discs
MIC	7240	Fluconazole 100 µg	5x50 Discs
MIC	7241	Fluconazole 25 µg	5x50 Discs
MIC	7242	Fluconazole 50 µg	5x50 Discs
MIC	7243	Flumequine 30 µg	5x50 Discs
MIC	7246	Framycetin 100 µg	5x50 Discs
MIC	7247	Furazolidone 100 µg	5x50 Discs
MIC	7248	Furazolidone 15 µg	5x50 Discs
MIC	7249	Furazolidone 50 µg	5x50 Discs
MIC	7250	Fusidic acid 30 µg	5x50 Discs
MIC	7251	Fusidic acid 5 µg	5x50 Discs
MIC	7252	Fusidic acid 50 µg	5x50 Discs
MIC	7254	Gatifloxacin 10 µg	5x50 Discs
MIC	7255	Gatifloxacin 30 µg	5x50 Discs
MIC	7256	Gentamicin 15 µg	5x50 Discs
MIC	7257	Gentamicin 200 µg	5x50 Discs
MIC	7258	Gentamicin 250 µg	5x50 Discs
MIC	7259	Gentamicin 50 µg	5x50 Discs
MIC	7260	Griseofulvin 10 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7262	Isepamicin 30 µg	5x50 Discs
MIC	7263	Isoniazid 1 µg	5x50 Discs
MIC	7264	Isoniazid 5 µg	5x50 Discs
MIC	7265	Itraconazole 50 µg	5x50 Discs
MIC	7266	Kanamycin 2 µg	5x50 Discs
MIC	7267	Kanamycin 5 µg	5x50 Discs
MIC	7268	Kanamycin 50 µg	5x50 Discs
MIC	7269	Ketoconazole 10 µg	5x50 Discs
MIC	7270	Levofloxacin 1 µg	5x50 Discs
MIC	7271	Levofloxacin 15 µg	5x50 Discs
MIC	7272	Lincomycin 10 µg	5x50 Discs
MIC	7273	Lincomycin 2 µg	5x50 Discs
MIC	7274	Lincomycin 30 µg	5x50 Discs
MIC	7275	Lincomycin / Neomycin 75 µg	5x50 Discs
MIC	7276	Lincomycin / Spectinomycin 109 µg	5x50 Discs
MIC	7277	Linezolid 5 µg	5x50 Discs
MIC	7278	Lomefloxacin 30 µg	5x50 Discs
MIC	7279	Marbofloxacin 5 µg	5x50 Discs
MIC	7280	Mecillinam 25 µg	5x50 Discs
MIC	7285	Metronidazole 2,5 µg	5x50 Discs
MIC	7286	Metronidazole 30 µg	5x50 Discs
MIC	7287	Metronidazole 4 µg	5x50 Discs
MIC	7288	Metronidazole 50 µg	5x50 Discs
MIC	7289	Metronidazole 80 µg	5x50 Discs
MIC	7290	Miconazole 50 µg	5x50 Discs
MIC	7291	Minocycline 25 µg	5x50 Discs
MIC	7292	Moxifloxacin 1 µg	5x50 Discs
MIC	7293	Mupirocin 20 µg	5x50 Discs
MIC	7294	Neomycin 5 µg	5x50 Discs
MIC	7295	Neomycin B 10 µg	5x50 Discs
MIC	7296	Netilmicin 10 µg	5x50 Discs
MIC	7297	Nitrocefin	5x50 Discs
MIC	7298	Nitrofurantoin 200 µg	5x50 Discs
MIC	7299	Nitrofurantoin 50 µg	5x50 Discs
MIC	7300	Nitroxoline 20 µg	5x50 Discs
MIC	7301	Nitroxoline 30 µg	5x50 Discs
MIC	7302	Norfloxacin 2 µg	5x50 Discs
MIC	7303	Norfloxacin 30 µg	5x50 Discs
MIC	7304	Novobiocin 10 µg	5x50 Discs
MIC	7305	Novobiocin 30 µg	5x50 Discs
MIC	7306	Penicillin / Novobiocin 40 µg	5x50 Discs
MIC	7309	Ofloxacin 2 µg	5x50 Discs
MIC	7310	Oleandomycin 15 µg	5x50 Discs
MIC	7312	Oxacillin 10 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7313	Oxacillin 5 µg	5x50 Discs
MIC	7314	Oxolinic acid 2 µg	5x50 Discs
MIC	7315	Oxytetracycline 30 µg	5x50 Discs
MIC	7317	Penicillin G 1.5 U µg	5x50 Discs
MIC	7318	Penicillin G 2 U µg	5x50 Discs
MIC	7319	Penicillin G 5 U µg	5x50 Discs
MIC	7322	Pipemidic acid 20 µg	5x50 Discs
MIC	7323	Pipemidic acid 30 µg	5x50 Discs
MIC	7324	Piperacillin 75 µg	5x50 Discs
MIC	7325	Piperacillin / Tazobactam 40 ug	5x50 Discs
MIC	7326	Piperacillin / Tazobactam 85 µg	5x50 Discs
MIC	7327	Pirlimycin 2 µg	5x50 Discs
MIC	7328	Polymixin B 100 U	5x50 Discs
MIC	7329	Polymixin B 50 U µg	5x50 Discs
MIC	7330	Posaconazole 5 µg	5x50 Discs
MIC	7331	Pristinamycin 15 µg	5x50 Discs
MIC	7332	PYR Test Discs	5x50 Discs
MIC	7333	PYR 50 Strips	50 strips
MIC	7334	Rifampin (Rifampicin) 15 µg	5x50 Discs
MIC	7335	Rifampin (Rifampicin) 2 µg	5x50 Discs
MIC	7336	Rifampin (Rifampicin) 30 µg	5x50 Discs
MIC	7339	Roxithromycin 15 µg	5x50 Discs
MIC	7340	Roxithromycin 30 µg	5x50 Discs
MIC	7341	Sarofloxacin 5 µg	5x50 Discs
MIC	7342	Sisomycin 10 µg	5x50 Discs
MIC	7343	Sodium Polyanethol Sulfonate	5x50 Discs
MIC	7344	Sparfloxacin 10 µg	5x50 Discs
MIC	7345	Spectinomycin 10 µg	5x50 Discs
MIC	7346	Spectinomycin 25 µg	5x50 Discs
MIC	7347	Spiramycin 100 µg	5x50 Discs
MIC	7348	Spiramycin 30 µg	5x50 Discs
MIC	7349	Streptomycin 25 µg	5x50 Discs
MIC	7350	Sulfadiazine 100 µg	5x50 Discs
MIC	7351	Sulfadiazine 300 µg	5x50 Discs
MIC	7352	Trimethoprim / Sulfadiazine 30 µg	5x50 Discs
MIC	7353	Sulfathiazole 300 µg	5x50 Discs
MIC	7354	Sulfisoxazole 30 ug	5x50 Discs
MIC	7355	Sulfisoxazole 300 µg	5x50 Discs
MIC	7356	Sulphamethoxazole 100 µg	5x50 Discs
MIC	7357	Sulphamethoxazole 200 µg	5x50 Discs
MIC	7358	Sulphamethoxazole 25 µg	5x50 Discs
MIC	7359	Sulphamethoxazole 300 µg	5x50 Discs
MIC	7360	Sulphamethoxazole 50 µg	5x50 Discs
MIC	7361	Terbinafine 30 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7362	Tetracycline 10 µg	5x50 Discs
MIC	7363	Thiamphenicol 30 ug	5x50 Discs
MIC	7364	Tilmicosin 15 µg	5x50 Discs
MIC	7365	Tobramycin 30 µg	5x50 Discs
MIC	7366	Trimethoprim 1,25 µg	5x50 Discs
MIC	7367	Trimethoprim 10 µg	5x50 Discs
MIC	7368	Trimethoprim 2,5 µg	5x50 Discs
MIC	7369	Trimethoprim 30 µg	5x50 Discs
MIC	7370	Trimethoprim 5 µg	5x50 Discs
MIC	7371	Triple Sulphonamides 1000 µg	5x50 Discs
MIC	7372	Triple sulphonamides 250 µg	5x50 Discs
MIC	7373	Triple sulphonamides 300 µg	5x50 Discs
MIC	7374	Tulathromycin 30 µg	5x50 Discs
MIC	7375	Tylosin 15 µg	5x50 Discs
MIC	7376	Tylosin 30 µg	5x50 Discs
MIC	7377	Tylosin 60 µg	5x50 Discs
MIC	7378	V Factor Disc	5x50 Discs
MIC	7379	Vancomycin 10 µg	5x50 Discs
MIC	7380	Vancomycin 5 µg	5x50 Discs
MIC	7381	Virginamycin 15 µg	5x50 Discs
MIC	7382	Voriconazole 1 µg	5x50 Discs
MIC	7383	Neomycin 10 µg	5x50 Discs
MIC	7384	Amoxicilin 20 ug	5x50 Discs
MIC	7385	Dicloxacilin 1 µg	5x50 Discs
MIC	7387	Blank Disc	5x50 Discs
MIC	7388	Meropenem / Cloxacillin 760 µg	5x50 Discs
MIC	7389	Meropenem / Dipicilonic acid 1010 µg	5x50 Discs
MIC	7390	Caspofungin 5 µg	5x50 Discs
MIC	7392	Cefdinir 50 µg	5x50 Discs
MIC	7393	Cefotaxime/Cloxacillin 780 µg	5x50 Discs
MIC	7394	Ceftazidime/Cloxacillin 780µg	5x50 Discs
MIC	7395	Ertapenem/Cloxacillin 760 µg	5x50 Discs
MIC	7396	Meropenem / Phenyl boronic acid 400 µg	5x50 Discs
MIC	7397	Meropenem 10 µg/boronic acid 610 µg	5x50 Discs
MIC	7398	Tiamullin 30 µg	5x50 Discs
MIC	7399	Trovafloxacin 10 µg	5x50 Discs
MIC	7400	Ampicillin 20 µg	5x50 Discs
MIC	7401	Carbenicillin 10 µg	5x50 Discs
MIC	7402	Cefepime 10 µg	5x50 Discs
MIC	7403	Cefepime 5 µg	5x50 Discs
MIC	7404	Cefepime/Clavulanic acid 15 µg	5x50 Discs
MIC	7405	Cefixime/Cloxacillin 40 µg	5x50 Discs
MIC	7406	Cefmetazole 30 µg	5x50 Discs
MIC	7407	Cefotaxime/Clavulanic Acid 15 µg	5x50 Discs

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7408	Cefotaxime/Sulbactam 40 µg	5x50 Discs
MIC	7409	Cefotaxime/Tazobactam 40 µg	5x50 Discs
MIC	7410	Cefoxitin / Cloxacillin 30 µg	5x50 Discs
MIC	7411	Ceftaroline 5 µg	5x50 Discs
MIC	7412	Ceftazidime/Clavulanic acid 20 µg	5x50 Discs
MIC	7413	Ceftriaxone/Tazobactam 40 µg	5x50 Discs
MIC	7414	Cephalothin 10 µg	5x50 Discs
MIC	7415	Clavulanic acid 15 µg	5x50 Discs
MIC	7416	Clotrimazole 20 µg	5x50 Discs
MIC	7417	Erythromycin 60 µg	5x50 Discs
MIC	7418	Fosfomycin 30 µg	5x50 Discs
MIC	7419	Kanamycin/cephalexin 45 µg	5x50 Discs
MIC	7420	Miconazole 30 µg	5x50 Discs
MIC	7421	Netilmicin 5 µg	5x50 Discs
MIC	7422	Oxolinic acid 10 µg	5x50 Discs
MIC	7423	Pipemidic acid 2 µg	5x50 Discs
MIC	7424	Rifampin (Rifampicin) 1 µg	5x50 Discs
MIC	7425	Tetracycline 5 µg	5x50 Discs
MIC	7426	Tiamullin 50 ug	5x50 Discs
MIC	7427	Azlocillin 30 µg	5x50 Discs
MIC	7429	Cefotiam 30 µg	5x50 Discs
MIC	7430	Colistin 40 µg	5x50 Discs
MIC	7433	Cephaloridine 30 µg	5x50 Discs
MIC	7434	Teicoplanin/norfloxacin	5x50 Discs
MIC	7435	Meropenem/EDTA 930 µg	5x50 Discs
MIC	7436	Ceftazidime / Avibactam 14 µg	5x50 Discs
MIC	7437	Cloxacillin 25µg	5x50 Discs
MIC	7438	Ertapenem/boranic acid 310µg	5x50 Discs
MIC	7500	Amikacin (0.016-256)	5x5 strips
MIC	7501	Amphotericin B antifungal (0.002-32)	5x5 strips
MIC	7502	Ampicillin (0.016-256)	5x5 strips
MIC	7503	Ampicillin/sulbactam (0.016-256)	5x5 strips
MIC	7504	Amoxicillin/clavulanic acid (0.016-256)	5x5 strips
MIC	7505	Azithromycin (0.016-256)	5x5 strips
MIC	7506	Aztreonam (0.016-256)	5x5 strips
MIC	7508	Cefaclor (0.016-256)	5x5 strips
MIC	7509	Cefazolin (0.016-256)	5x5 strips
MIC	7510	Cefepime (0.016-256)	5x5 strips
MIC	7511	Cefixime (0.016-256)	5x5 strips
MIC	7512	Cefoperazone/sulbactam (0.016-256)	5x5 strips
MIC	7513	Cefotaxime (0.002-32)	5x5 strips
MIC	7514	Cefotetan (0.016-256)	5x5 strips
MIC	7515	Cefoxitin (0.016-256)	5x5 strips
MIC	7517	Cefpodoxime (0.016-256)	5x5 strips

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7518	Ceftazidime (0.016-256)	5x5 strips
MIC	7519	Cefibuten (0.002-32)	5x5 strips
MIC	7520	Cefizoxime (0.002-32)	5x5 strips
MIC	7521	Ceftriaxone (0.002-32)	5x5 strips
MIC	7522	Cefuroxime (0.016-256)	5x5 strips
MIC	7523	Cephalothin (0.016-256)	5x5 strips
MIC	7524	Chloramphenicol (0.016-256)	5x5 strips
MIC	7525	Ciprofloxacin (0.002-32)	5x5 strips
MIC	7526	Clarithromycin (0.016-256)	5x5 strips
MIC	7527	Clindamycin (0.016-256)	5x5 strips
MIC	7528	Colistin (0.016-256)	5x5 strips
MIC	7529	Daptomycin (0.016-256)	5x5 strips
MIC	7530	Doripenem (0.002-32)	5x5 strips
MIC	7531	Doxycycline (0.016-256)	5x5 strips
MIC	7532	Enrofloxacin (0.002-32)	5x5 strips
MIC	7533	Ertapenem (0.002-32)	5x5 strips
MIC	7534	Erythromycin (0.016-256)	5x5 strips
MIC	7535	Fusidic acid (0.016-256)	5x5 strips
MIC	7536	Fluconazole antifungal (0.016-256)	5x5 strips
MIC	7537	Gatifloxacin (0.002-32)	5x5 strips
MIC	7538	Gentamicin (0.016-256)	5x5 strips
MIC	7540	Itraconazole antifungal (0.002-32)	5x5 strips
MIC	7541	Kanamycin (0.016-256)	5x5 strips
MIC	7542	Ketoconazole antifungal (0.002-32)	5x5 strips
MIC	7543	Levofloxacin (0.002-32)	5x5 strips
MIC	7544	Linezolid (0.016-256)	5x5 strips
MIC	7545	Meropenem (0.002-32)	5x5 strips
MIC	7546	Mecillinam (0.016-256)	5x5 strips
MIC	7547	Minocycline (0.016-256)	5x5 strips
MIC	7548	Moxifloxacin (0.002-32)	5x5 strips
MIC	7549	Nalidixic acid (0.016-256)	5x5 strips
MIC	7550	Netilmicin (0.016-256)	5x5 strips
MIC	7551	Nitrofurantoin (0.032-512)	5x5 strips
MIC	7552	Norfloxacin (0.016-256)	5x5 strips
MIC	7553	Ofloxacin (0.002-32)	5x5 strips
MIC	7554	Oxacillin (0.016-256)	5x5 strips
MIC	7555	Penicillin (0.016-256)	5x5 strips
MIC	7556	Piperacillin (0.016-256)	5x5 strips
MIC	7557	Piperacillin/tazobactam Antifungal (0.016-256)	5x5 strips
MIC	7558	Posaconazole (0.002-32)	5x5 strips
MIC	7559	Rifampicin (0.002-32)	5x5 strips
MIC	7560	Streptomycin (0.016-256)	5x5 strips
MIC	7561	Teicoplanin (0.016-256)	5x5 strips
MIC	7562	Tetracycline (0.016-256)	5x5 strips

CATEGORY	CAT.	PRODUCT	PACK SIZE
MIC	7563	Ticarcillin/clavulanic acid (0.016-256)	5x5 strips
MIC	7564	Tigecycline (0.016-256)	5x5 strips
MIC	7565	Tobramycin (0.016-256)	5x5 strips
MIC	7566	Trimethoprim (0.002-32)	5x5 strips
MIC	7567	Trimethoprim/sulfamethoxazole (0.002-32)	5x5 strips
MIC	7568	Vancomycin (0.016-256)	5x5 strips
MIC	7569	Fosfomicin (0.064 - 1024)	5x5 strips
MIC	7571	Cefotaxime (0.016 - 256)	5x5 strips
MIC	7572	Penicillin G (0.002-32)	5x5 strips
MIC	7573	Pefloxacin (0.002 - 32)	5x5 strips
MIC	7574	Streptomycin (0.064 - 1024)	5x5 strips
MIC	7575	Cefepime/Cefepime + Clavulanic Acid (0.25 -16 / 0.064 - 4)	5x5 strips
MIC	7576	Cefotaxime/Cefotaxime + Clavulanic acid (0.25 -16 / 0.016 -1)	5x5 strips
MIC	7577	Ceftazidime/Ceftazidime + Clavulanic Acid AMPC (0.5-32/0.064 - 4)	5x5 strips
MIC	7578	Cefotetan/Cefotetan + Cloxacillin (0.5 =32/0.5-32)	5x5 strips



Monitoring.

The control of environmental microbiological quality is essential in patient care spaces. Condalab offers the full product range required for air and surface control, as well as medical devices, so that the safety of the health personnel and patients is not compromised during their stay.

CATEGORY	CAT.	PRODUCT	PACK SIZE
Environmental	1024	Sabouraud Dextrose Agar	500 g
Environmental	1052	Macconkey Agar	500 g
Environmental	1058	Mueller Hinton Agar	500 g
Environmental	1102	Cetrimide Agar Base	500 g
Environmental	1400	Brain Heart Infusion Broth (BHI Broth)	500 g
Environmental	1423	MRSA Chromogenic Agar Base	500 g
Environmental	2059	Sabouraud Dextrose Agar + Chloramphenicol + Gentamicin	500 g
Environmental	2062	ESBL Chromogenic Agar	500 g
Environmental	2077	Vancomycin-Resistant Enterococcus (VRE) Chromogenic Agar	500 g
Medical devices	1068	Trypticasein Soy Agar (TSA)	500 g
Medical devices	1071	R2A Agar	500 g
Medical devices	1108	Blood Agar Base	500 g
Medical devices	1224	Trypticasein Soy Broth (TSB)	500 g
Medical devices	2043	Middlebrook Agar	500 g



Stains and Dyes.

Dyes and stains are mainly used in qualitative procedures for microorganism differentiation and identification based on their specific characteristics. This use is of vital importance, and it is usually the first step, when identifying a microorganism to establish a reliable clinical diagnosis in microbiology and pathological anatomy.

Find among our products the one that meets your needs, from classic stains in bacterial identification, to those used in specialized techniques for mycobacterial analysis, such as the Ziehl-Nielsen stain.

CAT.	PRODUCT	PACK SIZE
4600	Gram Stain KIT	4 x 250 ml
4603	Methylene Blue	1 x 1L
4610	Oxalate Crystal Violet	1 x 1L
4611	Chlorhidric Alcohol	1 x 1L
4612	Lugol	1 x 1L
4614	Acetone Ethanol Decolorant	1 x 1L
4616	Safranin	1 x 1L
4618	Fenicade Fuchsine	1 x 1L
4632	Thiazine Red	1 x 250 ml
4647	Auramide	1 x 1L
4651	Thiazine Red	1 x 1L
4655	Permanganate	1 x 1L
4687	CondaSafe Stain	1 ml
5005	Auramide	1 x 250 ml
5014	Permanganate	1 x 250 ml
5024	Chlorhidric Alcohol	1 x 250 ml
5058	Methylene Blue	1 x 250 ml
5059	Acetone Ethanol Decolorant	1 x 250 ml
5060	Fenicade Fuchsine	1 x 250 ml
5061	Lugol	1 x 250 ml
5062	Oxalate Crystal Violet	1 x 250 ml
5063	Safranin	1 x 250 ml



Other Reagents.

Since our founding in 1960, Condalab has become a specialist in the field of microbiology, which has positioned us today as a leading international privately held company.

However, Condalab is also recognized for supplying key ingredients such as agars, peptones and agaroses, among other products, to provide a complete and global solution according to our clients' needs.

CATEGORY	CAT.	PRODUCT	PACK SIZE
Agaroses	8001	Agarose MS - 6 Metagel	50, 100, 250 and 500 g
Agaroses	8010	Agarose D1 Low EEO	50, 100, 250 and 500 g
Agaroses	8017	Agarose D1 Low EEO GQT	50, 100, 250 and 500 g
Agaroses	8019	Agarose D1 Medium EEO	50, 100, 250 and 500 g
Agaroses	8024	Agarose D1 High EEO	50, 100, 250 and 500 g
Agaroses	8032	Agarose D2	50, 100, 250 and 500 g
Agaroses	8045	Agarose D5	50, 100, 250 and 500 g
Agaroses	8050	Agarose LM	50, 100, 250 and 500 g
Agaroses	8065	Agarose MS - 8	50, 100, 250 and 500 g
Agaroses	8067	Agarose MS - 12	50, 100, 250 and 500 g
Agaroses	8075	Agarose MS - 4	50, 100, 250 and 500 g
Agaroses	8090	Agarose FP DNA	50, 100, 250 and 500 g
Agaroses	8091	Agarose LM GQT	50, 100, 250 and 500 g
Agaroses	8092	Agarose LM Sieve	50, 100, 250 and 500 g
Agaroses	8093	Agarose NovaGel GQT	50, 100, 250 and 500 g
Agaroses	8100	Agarose E	50, 100, 250 and 500 g
Antibiotics	6805	Gentamicin	5 g
Antibiotics	6806	Kanamycin (Monosulphate)	5 g
Antibiotics	6807	Kanamycin (Monosulphate)	25 g
Strain Preservation	6082	Cryopreserve Beads Blue	50 cryotubes
Strain Preservation	6083	Cryopreserve Beads Green	50 cryotubes
Strain Preservation	6084	Cryopreserve Beads Yellow	50 cryotubes
Strain Preservation	6085	Cryopreserve Beads Natural	50 cryotubes
Strain Preservation	6086	Cryopreserve Beads Mixed	50 cryotubes
Strain Preservation	6087	Cryopreserve Beads Red	50 cryotubes
Ingredients	1612	Tryptone	500 g
Ingredients	1616	Bacteriological Peptone	500 g
Ingredients	1702	Yeast extract	500 g

CATEGORY	CAT.	PRODUCT	PACK SIZE
Ingredients	1800	European Bacteriological Agar	500 g
Ingredients	1802	American Bacteriological Agar	500 g
Ingredients	1806	Purified Agar	500 g
Ingredients	1900	Dextrose	500 g
Ingredients	1906	Sucrose	500 g
Ingredients	1935	Hemoglobin	500 g



Condalab

Inspired by knowledge

export@condalab.com | www.condalab.com

**For more information about the products and techniques
for clinical microbiology, do not hesitate to contact us.**