

INFRARED STERILIZER



Infrared Sterilizer



Sterilization only takes
5-7 seconds



Can be used in an anaerobic
environment



No open-flame, more safe



Easy to clean

TECHNICAL DATA

Model	DS-900	DS-900S
Power	150W	240W
Max. Temperature Inside The Heating Hole	825°C±50°C	
Heating Temperature	Two temperatures available 400°C/850°C	
Max. Diameter Of Sterilizer Goods (mm)	φ14	φ35
Length Of Heating Zone (mm)	150	100
Warm-up Time	10min	