

# VACUUM PUMP

---



## Oiless Diaphragm Vacuum Pump



Oil-free and no pollution



Illuminated switch



Overload,overheat and short circuit protection



High precision vacuum gauge

### TECHNICAL DATA

Model	VP-10L	VP-15L
Type	Negative pressure	
Suction Flow Rate	10L/min	15L/min
Limit Vacuum Degree	0.075Mpa	0.08Mpa
Noise	<40db	<45db
Enclosure Ratings	IP21	
Power Supply	20W	30W
Diaphragm	Rubber	
Protection Devices	Overload protection, over temperature protection, short circuit protection	
Gross/Net Weight	1.95/1.58 kg	

## Adjustable Speed Oil-free Diaphragm Vacuum Pump



The suction speed is adjustable



Illuminated switch



Overload, overheat and short circuit protection




High precision vacuum gauge

### TECHNICAL DATA

Model	VP-15S
Type	Negative pressure
Suction Flow Rate	0-15 L/min
Limit Vacuum Degree	0.08Mpa
Noise	< 45db
Enclosure Ratings	IP21
Power Supply	30W
Diaphragm	Rubber
Protection Devices	Overload protection, over temperature protection, short circuit protection
Gross/Net Weight	1.95/1.58 kg

## Oilless Diaphragm Vacuum Pump - 30L



Positive and negative pressure



Illuminated switch



Overload, overheat and short circuit protection



High precision vacuum gauge

### TECHNICAL DATA

Model	VP-30L
Type	Positive and negative pressure
Suction Flow Rate	30 L/min
Limit Vacuum Degree	0.085Mpa
Greatest Pressure	5.0bar
Noise	< 50db
Enclosure Ratings	IP21
Power Supply	50W
Diaphragm	Rubber
Protection Devices	Overload protection, over temperature protection, short circuit protection

## 6L Mini Vacuum Pump



9-speed speed regulation



Digital display speed control knob



Mini portable, strong power



Steady motion

### TECHNICAL DATA

Model	VP-6L	Pumping Rate	6L / Min
Type	Negative pressure suction	Ultimate Vacuum	0.06Mpa
Control Type	Digital display knob · 9 speeds	Power	10W

## Vacuum Filtration Apparatus



### TECHNICAL DATA

Receiving Flask	Sand Core Diameter	Sand Core Aperture
250ml	20mm	10um
500ml	40mm	10um
1000ml	40mm	10um
2000ml	40mm	10um
5000ml	40mm	30-50um
10000ml	40mm	30-50um



Boro 3.3



Tight sealing



Suitable for bottles with GL-45



Large diameter·Fast suction



### TECHNICAL DATA

Receiving Flask	Sand Core Diameter	Filter Membrane	Sand Core Aperture
2000ml	80mm	90-100mm	10um
5000ml	80mm	90-100mm	10um
10000ml	80mm	90-100mm	10um