

WATER BATH



Water Bath With Magnetic Stirrer



PID temperature control



Brushless Motor



Stainless steel inner tank



Strong Stirring Power

TECHNICAL DATA

Model	WB-6S	WBS-6Pro
Speed / Increment	/	100-1200RPM / 50RPM
Temperature Range	Room Temp. +5°C-100°C	Room Temp. +5°C-100°C
Liner Size	200*200*130mm(depth)-2L flask or 3L beaker	
Product Size (mm)	260*265*325	260*265*325
Display Accuracy	0.1°C	0.1°C
Power	800W	800W
Timing Range	1min-23h59min	1min-23h59min

Water Bath



Digital display thermostatic



High precise temperature control



One-piece stainless steel tank



High corrosion resistance

TECHNICAL DATA

Model	WB100-1 WB100-1F	WB100-2 WB100-2F	WB100-4 WB100-4F	WB100-6 WB100-6F	WB100-8F
Temperature Range	Room temp- 100°C				
Power	300W	600W	800W	1200W	1500W
Accuracy	±0.1°C				
Working Chamber Size (mm)	160*160*130	300*150*140	300*300*140	450*300*130	620*300*120

* WB100-F series are anti-dry series